

User and Care Instructions – Platypus Wheelchair



- Arm rest should only be lifted vertically, as shown in the above photo.
- Never stand on the footplates, this could cause the chair to tip and plates to break.
- Never stand on the frame of the chair. This will cause the chair to break.
- Never exceed the Safe Working Load. Every chair is supplied with stickers indicating the SWL. If this is removed see administration or call Para Mobility.
- Should the chair be damaged caused by misuse, it will not be covered under warranty
- When wheeling the chair, always go backwards over a step or bump, never step on the back of wheelchair to bring the front wheels up.

Cleaning your Platypus Pool / Wet Area Wheelchair

- Rinse your chair with fresh water after use to ensure all chlorine or salt water, grit, dirt and sand has been removed.
- When required, wash with soapy water. A detergent suitable for washing dishes will be suitable for cleaning all components of the chair, including the plastic frame, back rest and seat mesh, seat belt and lifting straps, wheels, brakes and any other component of the chair. Make sure to rinse all detergent off before using the chair in a pool.
- Wheels do not require any grease or oil.